



BASIC LESSONS
FOR UNDERSTANDING
EQUIPMENT



Raised Toilet Seat

Model #: _____

Provider Name: _____

Serial #: _____

Phone #: _____

Date: _____

Contact Person: _____

Setting up your equipment:

- Make sure the raised toilet seat fits your toilet. Some toilets may need modifications.
- Lift the toilet seat and place the raised toilet seat down. Inspect raised toilet seat for fit and stability prior to using it.
- The seat should rest flat on the bowl when properly installed.

Using your equipment:

- When using the raised toilet seat, you should center your weight as you sit. The seat may tip if you lean forward or side to side.
- It is recommended to use safety rails or a grab-bar for assistance.
- Always check to make sure that the raised toilet seat is secure and in place before use.

Maintenance:

- When cleaning the raised toilet, use a non-abrasive cloth and detergent with warm water. Clean the seat as needed.

Safety Issues:

- Observe the weight limit of your product.
- Users with limited physical capabilities should be supervised or assisted when using the seat.
- Always check that the raised toilet seat sits flat on the toilet bowl before it is used.

Please contact your equipment provider for any questions, adjustments, or repairs.

Your local provider:



©Copyright The VGM Group
B.L.U.E. Sheet images, drawings, and text
courtesy of VGM contracted vendors.

All Medicare DMEPOS suppliers must be in compliance with CMS Supplier Standards in order to obtain and retain their billing privileges. These standards, in their entirety, are listed in 42 C.F.R. pt. 424, sec 424.57(c). A supplier must disclose these standards to all customers/patients who are Medicare beneficiaries (standard 16). Supplier Standard 12 includes "A supplier ...must instruct beneficiaries on use of Medicare covered items".